

## ENSURE360 PULL/PUSH STATION

An integral part of the Sentrics Ensure360<sup>SM</sup> platform are the Push Station and Pull Station wall mounted emergency call stations. Push Station is a touch pad station offering two buttons for interaction, the face of the station and a button located on the side. Pull Station offers the same two buttons, and includes a pull cord hanging from the bottom. Both stations have a series of sensors, such as a motion detector to detect motion and ambient light to detect brightness. Other sensors include: temperature, barometric pressure, and humidity.

When a resident or caregiver pushes one of the buttons, the Push Station or Pull Station will communicate the desired action to the Extender network so that an emergency call or Assurance check can be completed.

### Installation Instructions

Without backplate (New install):

1. Remove Pull Station or Push Station from box
2. Twist to remove the back and separate the pieces
3. Measure 44" above finished floor (AFF), this is where the bottom of the back will rest
4. With the back placed properly, use the provided wall dogs to screw it into wall
5. Once mounted, remove the battery tab from the station and align the arrows on the front and back, twist until locked
6. Test station

### SENTRICS

245 Cedar Sage Drive  
Suite 250, Building 1B  
Garland, TX 75040  
[www.sentrics.net](http://www.sentrics.net)



### Specifications

**Dimensions:** 3.0" x 1.037"

## **ENSURE360 PULL/PUSH STATION**

With backplate (Retrofit):

1. Remove Pull Station or Push Station from box, get the backplate
2. Use the backplate to cover the existing hole from the previous station. The cord will hang about 42"
3. When the backplate is mounted properly, slide the station onto the backplate
4. Twist station to remove the front half, remove the battery tab from the station and align the arrows on the front and back, twist until locked
5. Test station

Parts:

- Pull.Carbon
- Push.Carbon
- 604106-Carbon retrofit backplate
- 301913 - 3pk CR2 batteries

### **Technical Specifications**

LED Patterns:

- Pull/Push powering on - Red, green, blue flash sequence



- Pull/Push transmitting an alarm - Blue flash
  - Pull/Push being updated - Green and yellow flash
- Supervision interval - 15 minutes  
Estimated battery life - Two years or longer

## **ENSURE360 PULL/PUSH STATION**

### **THE ENSURE360 SOLUTION**

Ensure360 is the perfect combination of Real-Time Location System (RTLS), environmental sensors, wearables, motion detectors, and traditional emergency call equipment designed to create a bubble of protection for your residents. The system, comprised of: a Headend server, a Desk console, a wireless mesh network made up of Extender access points, the Director network coordinator, and an array of Locator location beacons, will communicate with various end devices around the community to create a modular system to fit your needs! The end devices for Ensure360 include: the Pendant for tracking resident and staff location and providing a way to call for help from anywhere, the corded Pull Station, the touch pad Push Station, the Universal Transmitter with 2 dry contacts, and the Motion Sensor. Ensure360 will provide you with the most powerful emergency call system you've ever used.