

Bedside 360

The Bedside 360 monitoring device allows higher-acuity residents and their caregivers to feel safe while residents are in their beds. With both active and passive alerting options, caregivers are notified immediately on their Ensure360 Ecall system when the bedside monitoring device is activated.

The device is supervised, notifying staff if maintenance attention is required. The attached peripheral is likewise supervised, and if it becomes detached from the call station, staff receive a distinct alert notifying of the removal.





Accessories

- Bedside Wall-Mount Station with a 1/4" phono plug
- 7ft. Call Cord with a 1/4" phono plug



Mount vertically with the phono jack at the bottom of the device.











