



## CCSA-BD10-WTSC WIRELESS PRESSURE SENSITIVE PAD WITH RESTORAL

Pressure sensitive pads may be used to detect when a resident or patient leaves their bed or their wheelchair. Pads may be used in the bed or on the floor next to the bed. Restoral feature helps to facilitate the management and scheduling of alarms.



**Specifications** 

**Dimensions:** P-5-22 24" x 30"

P-6-22 9" x 16"

**Jacket:** 1/4" N/O

Reset: Momentary pushbutton N/O